

DynaMesh®-HIATUS implants are specially developed for prosthetic hiatoplasty and serve to achieve the long-lasting support and stabilisation of the diaphragm in the region of the oesophageal hiatus.

DynaMesh®-HIATUS

When selecting the mesh size, ensure sufficient overlap!

DynaMesh®-HIATUS	07 cm x 12 cm	REF 2611610712	BX = 1 piece
		REF 2611610713	BX = 3 pieces
	08 cm x 13 cm	REF 2611610813	BX = 1 piece
		REF 2611610814	BX = 3 pieces

VI014xx	DynaMesh®-HIATUS - Animation: Surgical Treatment of Hiatal Hernia in Laparoscopic Technique https://youtu.be/DHYSLgaM4uc	
---------	--	--

Use and Properties

Product	DynaMesh®-HIATUS
Field of application	hiatal hernia
Surgical access	laparoscopic
Surgical technique	hiatal hernia surgery with implant
Fixation	sutures / adhesives / tacks*
Shape stability	●
Atraumatic selvedges	●
Visible technology	●
Polymer (monofilament)	PVDF
Biocompatibility	●
Ageing resistance	●
Dynamometry	●
Tear propagation resistance	●
No scar plate formation	●
Classification (Klinge's classification [8])	1a



* Tacks may only be used if injury to the pericardium can be ruled out with certainty.

● Applies to all product sizes

Shape Stability Under Load

Conventional mesh structures deform under load. Constriction of the mesh in the region of the hiatus may reduce the distance between mesh implant and oesophagus, eventually causing mesh erosion.


DynaMesh®-HIATUS is based on a sophisticated textile design with rectangular pores, which even under load retain a high degree of shape stability.



High Effective Porosity


Mesh implants tend to shrink after incorporation in vivo.

DynaMesh®-HIATUS has a high effective porosity, which ensures that the mesh implant is thoroughly incorporated. During incorporation, the use of the proven and highly biocompatible PVDF polymer ensures that scarring is kept to a minimum. The good incorporation of the mesh implants combined with less scarring leads to minimisation of mesh shrinkage and permanently high flexibility of the incorporated implant.



Smooth, Warp-Knitted Selvedges

If the mesh does come into contact with the oesophagus in spite of all measures to prevent it, **DynaMesh®-HIATUS** has smooth selvedges that can minimise the danger of mesh erosion.



Distributed by: **DAHLHAUSEN®**
 P.J. Dahlhausen & Co. GmbH
 Alles Gute fürs Krankenhaus
 Emil-Hoffmann-Straße 53
 D-50996 Köln, Germany
 Tel.: +49 (0) 2236 - 39 13-0
 Fax: +49 (0) 2236 - 39 13-109
 www.dahlhausen.de
 info@dahlhausen.de