


DynaMesh®-SIS implants are designed as a midurethral sling for soft tissue reinforcement of the pelvic floor as part of the surgical treatment of stress urinary incontinence caused by a hypermobile urethra and/or intrinsic sphincter deficiency.

DynaMesh®-SIS


DynaMesh®-SIS	01 cm x 50 cm	REF 2621210156	BX = 1 piece
		REF 2623210156	BX = 3 pieces
DynaMesh®-SIS soft	01 cm x 50 cm	REF 2621410156	BX = 1 piece
		REF 2623410156	BX = 3 pieces
DynaMesh®-SIS visible	01 cm x 50 cm	REF 2621470156	BX = 1 piece
		REF 2623470156	BX = 3 pieces




Use and Properties

Product	DynaMesh®-SIS	DynaMesh®-SIS soft	DynaMesh®-SIS visible
Field of application	stress urinary incontinence (SUI)		
Surgical access	transvaginal		
Surgical technique	TVT - retropubic - inside-out / TOT - transobturator - inside-out - outside-in		
Fixation	none		
Atraumatic selvages		●	
Shape stability		●	
Defined elasticity		●	
Visible technology	●	●	●
Polymer (monofilament)		PVDF	
Biocompatibility		●	
Ageing resistance		●	
Dynamometry		●	
Tear propagation resistance		●	
Classification (Klinge's classification [8])		1a	

● Applies to all product sizes
● Does not apply



Retropubic
(inside-out / bottom-up)



Transobturator
(inside-out & outside-in)



DynaMesh®-SIS implants are positioned using the inside-out technique in case of a retropubic tape position, and using the outside-in or inside-out technique in case of a transobturator tape position.

DynaMesh®-SIS implants have a thread on both ends of the sling, which aids fixation to the surgical instrument.

Several reusable instruments are available separately to assist the positioning of DynaMesh®-SIS implants:



DynaMesh®-ISR01:
Instrument for retropubic positioning of DynaMesh®-SIS implants through transvaginal access using the inside-out technique.



Diameter: 5 - 7 cm

DynaMesh®-IST01/-IST02/-IST03:
Instrument set consisting of two instruments (right and left side) for transobturator positioning through transvaginal access using the inside-out or outside-in technique.



DynaMesh®-IVT01:
Instrument for transobturator positioning through transvaginal access using the outside-in technique.

VI040xx	DynaMesh®-SIS - Animation: Surgical Treatment of Stress Urinary Incontinence - SUI - TVT 8/4 https://youtu.be/OtxlWxz83FA	
---------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--

Distributed by: **DAHLHAUSEN®**
 P.J. Dahlhausen & Co. GmbH
 Alles Gute fürs Krankenhaus
 Emil-Hoffmann-Straße 53
 D-50996 Köln, Germany
 Tel.: +49 (0) 2236 - 39 13-0
 Fax: +49 (0) 2236 - 39 13-109
 www.dahlhausen.de
 info@dahlhausen.de