


DynaMesh®-PRS implants have been specially developed for pelvic floor reconstruction, in laparoscopic or open surgical technique, and serve to support and stabilise fascial structures, connective tissue and ligaments. The implants are used in the surgical treatment of the vaginal or cervical stump prolapse, as well as in the treatment of a concomitant cystocele and/or rectocele.

DynaMesh®-PRS

DynaMesh®-PRS soft	¹⁾ 02 cm x 16 cm	REF 2621350216	BX = 1 piece
DynaMesh®-PRS soft	²⁾ 03 cm x 23 cm	REF 2621350323	BX = 1 piece
DynaMesh®-PRS visible	³⁾ 03 cm x 23 cm	REF 2621750323	BX = 1 piece
DynaMesh®-PRS visible	⁴⁾ 3.3 cm x 24 cm	REF 2621750424 REF 2621750425	BX = 1 piece BX = 10 pieces
DynaMesh®-PRS visible	⁵⁾ 04 cm x 20 cm	REF 2621750420 REF 2621750421	BX = 1 piece BX = 10 pieces
DynaMesh®-PRS soft	⁶⁾ 05 cm x 27 cm	REF 2621350527	BX = 1 piece




VI046xx	DynaMesh®-PRS - Animation: Colposacropexy https://youtu.be/z4-9SVrGH4c	
VI048xx	DynaMesh®-PRS - Animation: Colposacropexy https://www.youtube.com/watch?v=sWlxsVbr7XY	
VI067xx	DynaMesh® MRI - Animation: MRI Reconstruction with DynaMesh®-PRS visible https://youtu.be/DHTxiYm3f2c	

Use and Properties


Product	DynaMesh®-PRS soft	DynaMesh®-PRS visible
	02 cm x 16 cm ⁽¹⁾ 03 cm x 23 cm ⁽²⁾ 05 cm x 27 cm ⁽⁶⁾	03 cm x 23 cm ⁽³⁾ 3.3 cm x 24 cm ⁽⁴⁾ 04 cm x 20 cm ⁽⁵⁾
Field of application	vaginal/cervical stump prolapse, concomitant cystocele/rectocele	
Surgical access	laparoscopic / open	
Surgical technique	colposacropepy / cervicosacropepy unilateral	
Fixation on vagina / cervix	sutures	
Fixation on sacrum	sutures / tacks	
Atraumatic selvedges	●	
Shape stability	●	
Defined elasticity	●	
Visible technology	● (1,2,6)	● (3,4,5)
Polymer (monofilament)	PVDF	
Biocompatibility	●	
Ageing resistance	●	
Dynamometry	●	
Tear propagation resistance	●	
Classification (Klinge's classification [8])	1a	

Application Examples:




Colpo-/cervicosacropepy

- unilateral



Colpo-/cervicosacropepy

- unilateral
- anterior/posterior mesh plasty
(for concomitant
cystocele/rectocele)



Colpo-/cervicosacropepy

- unilateral

● Applies to all product sizes
● Does not apply

Distributed by: **DAHLHAUSEN®**
 P.J. Dahlhausen & Co. GmbH
 Alles Gute fürs Krankenhaus
 Emil-Hoffmann-Straße 53
 D-50996 Köln, Germany
 Tel.: +49 (0) 2236 - 39 13-0
 Fax: +49 (0) 2236 - 39 13-109
 www.dahlhausen.de
 info@dahlhausen.de