



Female Pelvic Organ Prolapse
Vaginal/Cervical Stump Prolapse
or Uterine Prolapse



DynaMesh®-PR implants have been specially developed for pelvic floor reconstruction, in laparoscopic or open surgical technique, and serve to support and stabilise fascial structures, connective tissue and ligaments.

DynaMesh®-PR

DynaMesh®-PR soft	04 cm x 23 cm	REF 2621500423	BX = 1 piece
		REF 2623500423	BX = 3 pieces
DynaMesh®-PR visible	04 cm x 23 cm	REF 2621700423	BX = 1 piece




Use and Properties

Product	DynaMesh®-PR soft	DynaMesh®-PR visible
Field of application	vaginal/cervical stump or uterine prolapse, concomitant cystocele/rectocele	
Surgical access	laparoscopic / open	
Surgical technique	colposacropexy / cervicosacropexy / hysterovacropexy unilateral	
Fixation on vagina / cervix	sutures	
Fixation on sacrum	sutures / tacks	
Atraumatic selvedges		●
Shape stability		●
Defined elasticity		●
Visible technology	●	●
Polymer (monofilament)		PVDF
Biocompatibility		●
Ageing resistance		●
Dynamometry		●
Tear propagation resistance		●
Classification (Klinge's classification [8])		1 a

● Applies to all product sizes
● Does not apply


DynaMesh®-PR implants are used in the surgical treatment of the vaginal/cervical stump or uterine prolapse, as well as in the treatment of a concomitant cystocele/rectocele.

Application Examples:




Colpo-/cervicosacropexy

- unilateral
- fixation on vaginal/cervical stump




Colpo-/cervicosacropexy

- unilateral
- fixation on vaginal/cervical stump and anterior mesh plasty for concomitant cystocele




Colpo-/cervicosacropexy

- unilateral
- fixation on vaginal/cervical stump and posterior mesh plasty for concomitant rectocele



Hysterosacropexy

- unilateral
- posterior cervical fixation



Hysterosacropexy

- unilateral
- posterior cervical fixation and posterior mesh plasty for concomitant rectocele

VI086xx	DynaMesh®-PR - Animation: Colposacropexy https://youtu.be/jgTDDBVyFPE	
---------	--	---

Distributed by: **DAHLHAUSEN®**

P.J. Dahlhausen & Co. GmbH
Alles Gute fürs Krankenhaus
Emil-Hoffmann-Straße 53
D-50996 Köln, Germany
Tel.: +49 (0) 2236 - 39 13-0
Fax: +49 (0) 2236 - 39 13-109
www.dahlhausen.de
info@dahlhausen.de